

Installation Instructions


MODEL SM-007

1. **Remove sling seat from wheelchair** Most wheelchair seats are attached with 8 screws (4 on each side) and finish washers. Since the screws and washers will not be used to install this drop seat, make sure you retain them with the sling seat.

2. **Position seat and straps**

▲
Locate the **BACK** sticker on the bottom of the seat. Position the seat with this sticker at the back of the wheelchair.

Loop both front and back straps around the seat frame tubing and through cam buckles as illustrated in diagram A.



3. **Determine height and angle of seat**

Place the cushion to be used in the seat and adjust the length of the straps to locate the seating surface at desired height. Note that the front and back seat straps can be adjusted to different lengths to achieve an angled seating surface. Secure buckles.

4. **Secure Straps**


Wrap the four pre-installed hook and loop strap retainers around the straps to prevent accidental shifting. Diagram B

5. **Install Alarm Box**

Install 9 volt battery supplied into the alarm box.

Connect the cord from the seat sensor to the alarm box. The hook and loop strap allows alarm box to be hung from hand grip or strap can be opened to install on any portion of the frame. An upholstery clip is also included. Always route the cord away from moving parts.

Firmly press the switch to the ON position (blue LED light remains lit) and the alarm should sound. The alarm should silence automatically each time the user is seated in the wheelchair.



CAUTION

This device is not a substitute for proper supervision.
Check fall risk patients frequently.
Change battery when LED is flashing.

Note: Disconnect and remove alarm box prior to wheelchair washing.

Manufactured by:
210 Innovations, LLC
210 Leonard Drive
Groton, CT 06340
800-210-2298
Made in the U.S.A.